## INTENDED USE

The Triage<sup>®</sup> Cardiac Panel is a fluorescence immunoassay to be used with the Triage<sup>®</sup> Meter for the quantitative determination of Creatine Kinase MB, myoglobin and troponin I in EDTA whole blood and plasma specimens. The test is used as an aid in the diagnosis of myocardial infarction (injury).

## SUMMARY AND EXPLANATION OF THE TEST

The diagnosis of acute myocardial infarction (AMI) in a patient presenting with chest pain is difficult in many cases. The three major criteria outlined by the World Health Organization for differentiating chest pain associated with AMI from chest pain due to other non-cardiac reasons are: 1) patient history in addition to physical examination, 2) electrocardiographic data, and 3) changes in serum protein markers associated with myocardial infarction. At least two of these criteria must be fulfilled to appropriately diagnose an AMI.

Frequently, physical examination cannot differentiate AMI from other cardiac abnormalities. The electrocardiogram is useful in diagnosing AMI but is limited because it is diagnostic in only approximately 50% of AMI patients. Typically, Q wave formation and changes in the ST segment, elevation or depression, are indicative of AMI. However, the results of the electrocardiogram must be considered with the physical examination and clinical history of the patient. The electrocardiogram may be normal initially even though the patient is truly presenting with AMI.

Blood protein markers play an important role in the differential diagnosis of AMI when other indicators may be negative or questionable. Markers used in the diagnosis of myocardial infarction include: creatine kinase (CK), the MB isoenzyme of creatine kinase (CK-MB), lactate dehydrogenase, more specifically, LDH-1 and LDH-2, myoglobin, isoforms of the MM and MB isoenzymes of CK and the structural proteins of the troponin complex, i.e., troponin T and troponin I.

Following an AMI, the appearance of protein markers in the blood results from cellular necrosis initiated by an ischemic event. Those proteins that are present in the highest concentrations and those that are most soluble appear in the blood first, e.g., myoglobin. The structural and mitochondrial proteins of the myocytes appear later following infarction, e.g., CK-MB and proteins of the troponin complex, including troponin I.

Myoglobin is a cytoplasmic, soluble, heme protein present in muscle cells having a molecular weight of approximately 17,000 Daltons. Because of its relatively small size, high cellular concentration, and cytoplasmic location, myoglobin is released earlier than other cardiac markers following cellular necrosis or injury. Blood concentrations of myoglobin increase above the reference range within the first 2 hours following injury, reaching a peak between 6 and 8 hours after the onset of symptoms. Myoglobin returns to baseline or normal concentrations within 20-36 hours after tissue damage. Myoglobin is present in all types of muscle cells. Therefore, its appearance in blood is not necessarily associated with myocardial injury. Blood myoglobin concentrations may be elevated as a result of a variety of conditions that produce muscle damage. These include trauma, ischemia, surgery, exercise and a variety of degenerative muscular diseases. In this regard, myoglobin has its greatest value in the exclusion of myocardial infarction in the early hours following chest pain. Due to the rapid increase in blood myoglobin concentrations, followed by moderately sustained clearance, the utility of myoglobin is limited to the first 2-30 hours following tissue injury. Nevertheless, myoglobin is particularly useful when the clinical history of the patient is known.

Creatine Kinase MB (CK-MB) is an 82,000 Dalton cytosolic enzyme that is present in high concentrations in the myocardium. This isoenzyme of creatine kinase is frequently used in the diagnosis of acute myocardial infarction. Typically, CK-MB increases above normal within the first 4-8 hours following acute myocardial infarction, reaching maximum concentrations between 12 and 24 hours and returning to normal in approximately 3 days. CK-MB, like myoglobin, is not specifically localized in cardiac muscle. Blood concentrations of CK-MB can be elevated as a result of acute or chronic muscle damage, including strenuous exercise and trauma. Nonetheless, measurements of blood CK-MB concentrations are widely relied on for the management of patients having an AMI.

The contractile proteins of the myofibril have gained increased popularity as cardiac specific markers for acute myocardial infarction and myocardial damage. These include two specific proteins of the contractile regulatory

complex, troponin I and troponin T. Troponin I and troponin T isolated from cardiac muscle have unique amino acid sequences that enable the development of specific antibodies to the cardiac proteins.

The amino terminal amino acid sequence of the cardiac isotype of troponin I has 31 amino acid residues that are not present in either of the two isoforms of troponin I in skeletal muscle. Therefore, immunoassays specific for cardiac troponin I are used in the evaluation of patients suspected of experiencing an AMI. Blood troponin I concentrations become elevated between 4 and 8 hours following an AMI. The concentration peaks between 12 and 16 hours and remains elevated for 5-9 days following damage to the myocardium. Cardiac troponin I is primarily elevated as a result of myocardial infarction. However, cardiac troponin I also may be elevated as a result of minor cardiac injury that includes: unstable angina, cardiac contusions, cardiac transplant, coronary artery bypass graft surgery, physical trauma to the heart, congestive heart failure and other conditions that may damage the myocardium. Moreover, cardiac troponin I does not appear to be elevated as a result of skeletal muscle injury. Due to the increased analytical specificity and the increased duration of elevation, cardiac troponin I has become an important marker in the diagnosis and evaluation of patients suspected of having an AMI. Simultaneous quantification of myoglobin, CK-MB and cardiac troponin I following AMI can greatly assist the physician in the management of patients suspected of presenting with an AMI.

Troponin I concentrations also have been described in the scientific literature to provide prognostic information related to the risk of future cardiac events and mortality in patients with acute coronary syndromes. More recently, it has been demonstrated that a multimarker analysis including troponin I, CK-MB, and myoglobin provides better risk stratification than a single-marker approach.

## **PRINCIPLES OF PROCEDURE**

The Triage<sup>®</sup> Cardiac Panel is a fluorescence immunoassay for the quantitative determination of the cardiac proteins, CK-MB, myoglobin and troponin I in whole blood and plasma specimens using EDTA as the anticoagulant. After addition of the sample to the sample port, the cells are separated from the plasma via a filter contained in the device. A predetermined quantity of plasma is allowed to react with fluorescent antibody conjugates within the reaction chamber. After sufficient incubation has occurred, the reaction mixture flows down the device

detection lane. Complexes of the analytes and fluorescent antibody conjugates are captured on discrete zones resulting in binding assays that are specific for each analyte. The concentration of the analyte in the specimen is directly proportional to the fluorescence detected.

## REAGENTS

The Triage<sup>®</sup> Cardiac Panel contains all the reagents necessary for the simultaneous quantification of the cardiac proteins CK-MB, myoglobin and troponin I in plasma and whole blood.

The device contains:

- Murine monoclonal and polyclonal antibodies against CK-MB, murine monoclonal and polyclonal antibodies against myoglobin, and murine monoclonal and goat polyclonal antibodies against troponin I labeled with a fluorescent dye and immobilized on the solid phase, and stabilizers.

## WARNINGS AND PRECAUTIONS

- For In Vitro Diagnostic Use.
- Carefully follow the instructions and procedures described in this insert.
- Patient specimens, used test devices and transfer pipettes may be potentially infectious. Proper handling and disposal methods should be established by the laboratory director in accordance with local, state and federal regulations.
- Proper laboratory safety techniques should be followed at all times when working with blood specimens because they are potentially infectious.
- The Triage<sup>®</sup> Cardiac Panel, like all cardiac marker tests, does not provide a definitive diagnosis of myocardial infarction. As with all in vitro diagnostic tests, the test results should be interpreted by the physician in conjunction with other laboratory test results and patient clinical findings.
- The transfer pipette should be used for one specimen only. Discard after single use.
- Keep the device in the sealed pouch until ready to use. Discard after single use.

## STORAGE AND HANDLING REQUIREMENTS

The Triage<sup>®</sup> Cardiac Panel is stable (while in the sealed pouch) up to the expiration date when stored refrigerated at 2° - 8°C (35° - 46°F). Once removed from refrigeration, the pouched Triage<sup>®</sup> Cardiac Panel is stable for

up to 14 days, but not beyond the expiration date printed on the pouch. If not used on the same day of removal from refrigeration, gently write the date of removal from the refrigerator and the date to discard on the foil pouch and/or the kit box (use a soft, felt tip marker).

- Once removed from refrigeration, allow a minimum of 15 minutes for device to reach room temperature while in the sealed pouch.
- Do not remove the device from the pouch until ready for use.

### **SAMPLE COLLECTION AND PREPARATION**

Whole blood or plasma using EDTA as the anticoagulant is the specimen required for testing with this product. Other blood specimen types have not been evaluated.

- Test blood specimens on the Triage® Cardiac Panel immediately or within 4 hours of sample collection. If testing cannot be completed within 4 hours, the plasma should be separated and stored at -20 °C until it can be tested.
- Transport specimens at room temperature or chilled and avoid extreme temperatures.
- The specimen must be at room temperature and be homogeneous before testing.
- Mix the blood specimen by gently inverting the tube several times before transferring blood to the test device.
- It is recommended to avoid using severely hemolyzed specimens whenever possible. If a specimen appears to be severely hemolyzed, another specimen should be obtained and tested.

### **MATERIALS PROVIDED**

Biosite® Triage® Cardiac Panel Catalog No. 97000HS

Kit contains:

Test Devices	25
Transfer Pipettes	25
Reagent Code Chip®	1
Printer Paper	1 roll

### **MATERIALS REQUIRED BUT NOT PROVIDED**

Triage® Meter

Triage® Cardiac Controls

## PROCEDURE

### USING THE TRANSFER PIPETTE

1. Depress the larger (top) bulb completely and insert the tip into the specimen.
2. Release the bulb slowly. The tube should fill completely. Some fluid must remain in the smaller (lower) bulb.
3. Place the tip of the pipette into the sample port of the device and depress larger bulb completely. The entire amount of sample in the pipette tube must be dispensed. The sample in the smaller (lower) bulb will not be expelled.
4. Remove the tip from the sample port and then release the bulb.

### TEST PROCEDURE

#### STEP 1. Add Sample to Triage® Cardiac Panel Device

Open the pouch, write the patient ID number on the label of the device, and add EDTA-anticoagulated whole blood or plasma into the sample chamber using the transfer pipette provided (refer to pipette instructions).

#### STEP 2. Insert Device into Triage® Meter

Allowing a moment for the sample to absorb into the sample port of the device, insert the Triage® Cardiac Panel device into the Triage® Meter. This step must be completed no longer than 30 minutes after step 1.

#### STEP 3. Read the Results on the Triage® Meter Display Screen

The test results are displayed when the test is complete. After the device has been analyzed, it is released by the Triage® Meter.

## RESULTS

The results of the Triage® Cardiac Panel are calculated automatically by the Triage® Meter. The concentrations of all three analytes are displayed on the screen. The operator has the option to print the results. For additional information, refer to the User Guide for the Triage® Meter.

## STANDARDIZATION

The Triage® Cardiac Panel has been standardized using purified protein preparations of CK-MB, myoglobin and cardiac troponin I based on the mass (concentration) of analyte present in EDTA-anticoagulated plasma.

## QUALITY CONTROL AND ACCEPTABILITY OF RESULTS

The use of positive and/or negative controls is recommended to assess each shipment of product, and at least every 30 days, or whenever the laboratorian wishes to verify the performance of the controls or the test. Biosite provides two external control solutions for this purpose. One contains analyte concentrations near the diagnostic cut-off. The other contains concentrations of analytes that are observed in samples from patients experiencing an AMI. It is recommended that these controls be used to test each new shipment or lot of product and to verify the system integrity and operator proficiency.

Each Triage® Cardiac Panel contains two internal controls that satisfy routine quality control requirements. These controls indicate that sufficient sample was applied to the panel, the unbound fluorescent label was washed sufficiently from the detection zone, and the panel was inserted and read properly by the meter. An unacceptable result from either control causes the display of a warning message on the Triage® Meter indicating that the test should be repeated.

For further information regarding the complete quality control of the system refer to the User Guide for the Triage® Meter.

The QC Device (may also be named QC Simulator) provided with the Triage® Meter should be tested each day of patient testing to verify instrument performance. Alternately, the QC Device should be tested upon set-up of the Triage® Meter and whenever required by local, state and federal regulations. Refer to the User Guide for a description of the QC Device.

## LIMITATIONS OF THE PROCEDURE

The results of the Triage® Cardiac Panel should be evaluated in the context of all the clinical and laboratory data available. In those instances where the laboratory results do not agree with the clinical evaluation, additional tests should be performed accordingly.

This test has been evaluated with whole blood and plasma using EDTA as the anticoagulant. Serum and blood or plasma specimens obtained using other anticoagulants have not been evaluated and should not be used (e.g. heparin or citrate).

There is the possibility that factors such as technical or procedural errors, as well as additional substances in blood specimens that are not listed below, may interfere with the test and cause erroneous results.

## PERFORMANCE CHARACTERISTICS

### ANALYTICAL SENSITIVITY

The analytical sensitivity or lowest detectable concentration that is distinguishable from zero for the three analytes was determined by testing a zero calibrator 20 times each using 3 lots of reagents and 5 meters on 3 days. The analytical sensitivity of each assay on the Triage® Cardiac Panel is presented below:

Troponin I:	0.05 ng/mL
CK-MB:	1.0 ng/mL
Myoglobin:	5 ng/mL

### MEASURABLE RANGES

Troponin I:	0.05 - 30 ng/mL
CK-MB:	1.0 - 80 ng/mL
Myoglobin:	5 - 500 ng/mL

## INTERFERING SUBSTANCES

Hemoglobin (up to 1,000 mg/dL), lipids (cholesterol up to 1,000 mg/dL and triglycerides up to 1,000 mg/dL) or bilirubin (up to 20 mg/dL) added to EDTA-anticoagulated plasma containing the three analytes did not interfere with the recovery of the analytes. These substances failed to produce a positive response in a sample that did not contain any of the analytes of interest. However, severely hemolysed specimens should be avoided whenever possible. When a sample appears to be severely hemolyzed, another specimen should be obtained and tested.

The hematocrit was varied between 30% and 60% and with no significant effect on the recovery of troponin I, CK-MB or myoglobin.

### PHARMACEUTICALS

The following drugs were evaluated for potential cross-reactivity and interference in the Triage® Cardiac Panel. All drugs were tested at concentrations that represent the blood concentrations that would result from a maximal therapeutic dose and at least twice the maximal therapeutic dose. None of the drugs interfered with the recovery of CK-MB, myoglobin or troponin I. Additionally, these drugs did not produce a significant response when tested in a specimen containing none of the analytes of interest. There was no significant interference with the analyte, nor was there any assay cross-reactivity.

Acetaminophen	Dopa, 1-alpha-methyl	Nitrofurantoin
Acetylsalicylic acid	Dopamine	Nitroglycerin
Allopurinol	Enalapril maleate	Oxazepam
Amiodarone	Erythromycin	Oxytetracyclin
Ampicillin	Furosemide	Phenobarbital
Ascorbic acid	Heparin	Phenytoin
Atenolol	Hydrochlorothiazide	Probenecid
Caffeine	Indomethacin	Procainamide
Captopril	Isosorbide dinitrate	Propranolol
Chloramphenicol	Lisinopril	Quinidine
Cyclosporine	Lovastatin	Sulfamethoxazole
Diclofenac	L-thyroxine	Theophylline
Digoxin	Nicotine	Trimethoprim
Diltiazem	Nicotinic acid	Verapamil
Dipyridamole	Niphedipine	Warfarin

## PROTEINS

Protein	ng/mL	Reactivity with Related Protein		
		CK-MB % Cross-reactivity%	Troponin I Cross-reactivity %	Myoglobin % Cross-reactivity
<b>Control</b>				
Actin	500	0.0%	0.0%	0.0%
Actin	1000	0.0%	0.0%	0.0%
CK-BB	15.6	0.3%	0.0%	0.0%
CK-BB	31.2	0.8%	0.0%	0.0%
CK-BB	62.5	0.9%	0.0%	0.0%
CK-BB	125	1.5%	0.0%	0.0%
CK-BB	250	2.4%	0.0%	0.0%
CK-BB	500	3.3%	0.0%	0.0%
CK-MM	250	0.2%	0.0%	0.0%
CK-MM	500	0.0%	0.0%	0.0%
CK-MM	5000	0.0%	0.0%	0.0%
cTnC	2000	0.0%	0.0%	0.0%
cTnT	2000	0.0%	0.0%	0.0%
Myosin	2000	0.0%	0.0%	0.0%
sTnI	500	0.0%	0.0%	0.0%
sTnI	1000	0.0%	0.0%	0.0%
sTnT	500	0.0%	0.0%	0.0%
sTnT	1000	0.0%	0.0%	0.0%
Tropomyosin	2000	0.0%	0.0%	0.0%

In addition to the proteins tested above, the Triage® Cardiac Panel was tested with regard to the ability of the troponin I test to detect various complexes of cardiac troponin I. The results below demonstrate that the Triage® Cardiac Panel recognizes 5 forms of cardiac troponin I on an equimolar basis.

#### Reactivity with Various Forms of Cardiac Troponin I

Troponin Form	Troponin Recovery (ng/mL)	Troponin Recovery (%)
Troponin I, Oxidized	1.21	100
Troponin I, Reduced	1.12	93
Troponin I-C Complex	1.52	125
Troponin I-T Complex	1.39	115
Troponin C-T-I Complex	1.19	99

Recent reports show that cardiac troponin I is released as binary and ternary complexes, in addition to free troponin I from patients suffering from AMI.

In light of these reports it would seem that assays for cardiac troponin I should be able to detect the analyte in each of its forms on an equimolar basis (free and complex).

## IMPRECISION

Within-day and total imprecision were determined using the ANOVA model by testing control materials and human plasma pools that had the respective analytes added at concentrations near the decision points of the assay and throughout the range of the standard curve. The study was conducted over 10 days, testing each control 10 times per day.

**CK-MB****Average Within Day Precision**

Mean (ng/mL)	SD (ng/mL)	CV
4.8	0.5	11.4%
15.8	2.1	13.4%
38.4	5.5	14.3%

**Average Total Precision**

Mean	SD	CV
4.8	0.6	11.6%
15.8	2.2	14.2%
38.4	5.4	14.1%

Data reflects 10 measurements per day for 10 days, with each device read on 1 meter

**MYOGLOBIN****Average Within Day Precision**

Mean (ng/mL)	SD (ng/mL)	CV
77.4	7.8	10.1%
111.8	11.6	10.4%
217.6	23.5	10.8%

**Average Total Precision**

Mean	SD	CV
77.4	9.0	11.6%
111.8	13.7	12.2%
217.6	28.2	13.0%

Data reflects 10 measurements per day for 10 days, with each device read on 1 meter

**TROPONIN I****Average Within Day Precision**

Mean (ng/mL)	SD (ng/mL)	CV
0.12	0.02	16.3%
0.22	0.03	11.7%
0.39	0.05	12.0%
1.8	0.2	11.5%
15.2	1.6	10.3%

**Average Total Precision**

Mean	SD	CV
0.4	0.05	12.0%
1.8	0.2	12.8%
15.2	1.7	11.2%

Data reflects 10 measurements per day for 10 days, with each device read on 1 meter

## EXPECTED VALUES

**HEALTHY VOLUNTEERS**

CK-MB and myoglobin concentrations were determined using specimens obtained from 452 apparently healthy individuals (264 women and 188 men) The 95th percentiles of concentrations for each analyte are shown below.

Analyte	95th Percentile
CK-MB	< 4.3 ng/mL
Myoglobin	<107 ng/mL

Troponin I concentrations were determined using specimens obtained from 323 apparently healthy individuals (168 women and 155 men). The 95th, 97.5th and 99th percentiles are shown below.

Analyte	95th Percentile	97.5th Percentile	99th Percentile
Troponin I	< 0.05ng/mL	<0.05 ng/mL	<0.05 ng/mL

## PATIENTS WITH SKELETAL MUSCLE INJURY AND RENAL DISEASE

Two additional groups of patients were evaluated for the presence of the various analytes. Both myoglobin and CK-MB are known to be potentially elevated in these conditions and diseases. Most of the patients' analyte concentrations were evaluated at a single time. Those patients who were suffering from renal insufficiency were evaluated at the time they received their dialysis. Cardiac involvement was not considered during the initial screening of the patients in the study. Some patients were determined to have cardiac injury or contusions after the initial diagnosis. Twelve of the twenty-one specimens having elevated troponin I concentrations obtained from patients having a primary diagnosis of skeletal muscle trauma were elevated as a result of cardiac involvement.

### SKELETAL MUSCLE INJURY

	CK-MB (ng/mL)	Troponin-I (ng/mL)	Myoglobin (ng/mL)
Number of patients/samples	117/189	117/189	117/189
Cut-Off	4.3	0.4	107
No. Samples Above Cut-Off	121	21	165
No. Samples from Patients with Cardiac Involvement	15	12	15
Clinical Specificity	83/189 x 100 44%	180/189 x 100 95%	39/189 x 100 21%

## RENAL PATIENTS

	CK-MB (ng/mL)	Troponin-I (ng/mL)	Myoglobin (ng/mL)
Number of patients	80	80	80
Cut-Off	4.3	0.4	107
No. Samples Above Cut-Off	22	5	74
No. Samples from Patients with Cardiac Involvement	5	5	5
Clinical Specificity	63/80 x 100 79%	80/80 x 100 100%	11/80 x 100 16%

Those conditions that result in myocardial cell damage potentially can cause increased blood concentrations of any of these analytes. For example, troponin I concentrations have been reported to be elevated as a result of unstable angina, congestive heart failure, myocarditis, and cardiac surgery, including invasive cardiac testing and cardiac contusions. Additionally, both CK-MB and myoglobin have been reported to be elevated in both skeletal muscle injury and renal disease.

## INTERPRETATION OF RESULTS

Temporal elevations of CK-MB, myoglobin and troponin I are observed in patients diagnosed with myocardial infarction. However, CK-MB and myoglobin, but not cardiac troponin I, may be elevated in renal disease and skeletal muscle injury. Cardiac troponin I appears to be elevated only in those diseases that directly involve the heart. Collectively, the diagnosis of myocardial infarction should include measurement of these cardiac related proteins and other clinical information including patient history and electrocardiographic data. Other conditions that may result in elevated cardiac proteins are: cardiac contusions, myocarditis, invasive examination of the heart, coronary artery bypass surgery, congestive heart failure and unstable angina. Therefore, these data must be considered when interpreting the results.

These values are representative. Each laboratory should establish a reference range that is representative of the patient population to be evaluated. Additionally, each laboratory should consider the current practice in the evaluation of patients experiencing chest pain and AMI at their respective institution.

## **CLINICAL PERFORMANCE IN THE EVALUATION OF CHEST PAIN**

The Triage® Cardiac Panel was evaluated in studies at 4 sites. In addition to the ranges of expected concentrations in apparently healthy individuals, patients with renal disease and patients suffering from acute muscle injury, the clinical sites evaluated patients for whom a diagnosis of myocardial infarction was indicated. Clinical diagnosis of myocardial infarction was based on satisfying at least two of the three criteria outlined below:

- Chest pain (discomfort) for a duration of at least 20 minutes
- Electrocardiographic changes that are consistent with myocardial infarction
- Temporal changes in cardiac enzymes (markers)

Those patients who did not satisfy 2 of the 3 criteria listed above were classified in the rule-out group.

The diagnostic sensitivity and specificity were evaluated by comparing the marker concentration to the discharge diagnosis for each patient. Inasmuch as the WHO criteria used in the diagnosis of myocardial infarction does not encompass the diagnosis of minor myocardial injury, the diagnostic specificity of troponin I may appear to be less than CKMB when using these criteria.

## **CLINICAL SENSITIVITY AND SPECIFICITY BY TIME INTERVAL**

Temporal elevations of all three cardiac markers are useful in the management of patients with chest pain, to aid in the diagnosis myocardial infarction and assessment of patients exhibiting chest pain. Serial sampling of the patient's blood following myocardial infarction is recommended, as changes in the marker concentrations may also be diagnostic. In particular, it has been demonstrated that temporal changes in myoglobin concentration provide additional diagnostic information that would otherwise not be identified if using a single time point. It is recommended that each hospital establish a suitable sampling protocol in addition to establishing an appropriate reference range. The cut-off concentrations for troponin I

(0.4 ng/mL), CK-MB (4.3 ng/mL) and myoglobin (107 ng/mL) were used to calculate the clinical sensitivity and specificity.



225 patients experiencing symptoms of acute myocardial infarction were evaluated. 207 specimens were obtained and evaluated from 72 patients diagnosed with myocardial infarction. An additional 316 samples from 153 patients were obtained and evaluated from patients for whom a diagnosis of AMI was excluded. Included were patients with unstable angina, coronary artery disease and other causes of chest pain but where AMI was ruled out.


# of samples	Clinical Sensitivity			
	0-6 hrs.	6-12 hrs.	12-24 hrs.	> 24 hrs.
Cardiac Troponin I Sensitivity	40	32	44	91
95% Confidence Interval	65%	72%	93%	96%
CK-MB Sensitivity	50.2%-79.8%	56.3%-87.5%	85.7%-100%	91.4%-99.8%
95% Confidence Interval	78%	78%	80%	85%
Myoglobin Sensitivity	64.6%-90.4%	63.8%-92.4%	67.6%-91.5%	77.2%-92.0%
95% Confidence Interval	75%	75%	73%	74%
	61.1%-88.4%	60.0%-90.0%	59.6%-85.9%	64.6%-82.7%

# of samples	Clinical Specificity			
	0-6 hrs.	Time		
		6-12 hrs.	12-24 hrs.	> 24 hrs.
	89	66	90	71
Cardiac Troponin I Specificity	100%	97.00%	94.40%	90.10%
95% Confidence Interval	100%-100%	92.8%-100%	89.7%-99.2%	83.2%-97.1%
CK-MB Specificity	91.00%	86.40%	82.20%	88.70%
95% Confidence Interval	85.1%-97.0%	78.1%-94.6%	74.3%-90.1%	81.4%-96.1%
Myoglobin Specificity	74.20%	81.80%	67.80%	71.80%
95% Confidence Interval	65.1%-83.3%	72.5%-91.1%	58.1%-77.4%	61.4%-82.3%

## ROC ANALYSIS OF TROPONIN I, CK-MB AND MYOGLOBIN

The graph below depicts the clinical sensitivity and specificity of cardiac troponin I, CK-MB and myoglobin when using various cut-off concentrations. The upper end of normal values were used as the cut-off for CKMB (4.3 ng/mL), myoglobin (107 ng/mL), and troponin I (0.4 ng/mL). Additionally these values were used as the cut-off concentrations for the statistics provided above. Each laboratory should establish their own diagnostic cut-off concentrations based on the clinical practice at their respective institutions.

Biosite's express and implied warranties (including implied warranties of merchantability and fitness) are conditioned upon observance of Biosite's published directions with respect to the use of Biosite's products.

UNDER NO CIRCUMSTANCES WHATSOEVER SHALL BIOSITE INCORPORATED BE LIABLE FOR ANY INDIRECT OR CONSEQUENTIAL DAMAGES.

In the U.S. and Canada, for Technical Assistance call  
Biosite® Technical Services at  
1-888-BIOSITE (246-7483)

Manufacturer:  
Biosite Incorporated  
11030 Roselle St.  
San Diego, California 92121  
+1 (858) 455-4808  
[www.biosite.com](http://www.biosite.com)

Made in USA

Purchase of this product licenses its use under U.S. Patent Numbers 5458852, 5763189, 5885527, 6019944, 6057098, 6074616, 6143576, 6156270, 6194222, 6238931, and 6251687. Additional patents pending.

Triage® and Biosite® are registered trademarks of Biosite Incorporated.

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AHA Medical/Scientific Statement, ACC/AHA Guidelines for the Early Management of Patients with Acute Myocardial Infarction. *Circulation* **82**: 664-707, 1990.

Bodor, G.S., Porter S., Landt, Y. and Ladenson, J.H. Development of Monoclonal Antibodies Specific for Troponin I and Preliminary Results in Suspected Cases of Myocardial Infarction. *Clin. Chem.* **38**: 2203-2214, 1992.

Puleo, P.R., Guadagno P.A., Roberts, R., Scheel, M.V., Marian, A.J., Churchill, D., and Perryman, M.B. Early Diagnosis of Myocardial Infarction for Subforms of Creatine Kinase-MB. *Circulation* **82**: 759-764, 1990.

Marin, M.M., and Teichman, S.L. Use of Rapid Serial Sampling of Creatine Kinase MB for Very Early Detection of Myocardial Infarction in Patients with Acute Chest Pain. *Am. Heart J.* **123**: 354-361, 1992.

Gerhardt, W., Waldenstrom, J., Horder, M., Hofvendahl, S., Billstrom, R., Ljungdahl, R., Berning, H., and Bagger, P. Creatine Kinase and Creatine Kinase B-Subunit Activity in Serum in Cases of Suspected Myocardial Infarction. *Clin. Chem.* **26**: 277-283, 1982.

Lee, T.H. and Goldman, L. Serum Enzyme Assays in the Diagnosis of Acute Myocardial Infarction: Recommendations Based on Quantitative Analysis. *Ann. Int. Med.* **105**: 221-233, 1986.

Vaidya, H.C., Maynard, Y., Dietzler, D.N., and Ladenson, J.H. Direct Measurement of Creatine Kinase-MB Activity in Serum after Extraction with a Monoclonal Antibody Specific to the MB isoenzyme. *Clin. Chem.* **32**: 657-663, 1986.

Hedges, J.R., Rouan, G.W., Tolzis, R., Goldstein-Wayne, B., and Stein, E.A. Use of Cardiac Enzymes Identifies Patients with Acute Myocardial Infarction Otherwise Unrecognized in the Emergency Department. *Ann. Emerg. Med.* **16**: 248-252, 1987.

Apple, F.S. Diagnostic Use of CK-MM and CK-MB Isoforms for Detecting Myocardial Infarction. *Clin. Lab. Med.* 9: 643-655, 1989.

Hedges, J.R., Swanson, J.R., and Heeter, C. Prospective Assessment of Presenting Serum Markers for Cardiac Risk Stratification. *Ac. Emerg. Med.* 3: 27-33, 1996.

Willerson, J.T., Clinical Diagnosis of Acute Myocardial Infarction. *Hosp. Prac.* 24: 65-77, 1989.

Cummins, B., Auckland, M.S., and Cummins, P. Cardiac-Specific Troponin I Radioimmunoassay in the Diagnosis of Acute Myocardial Infarction. *Am. Heart J.* 113: 1333-1344, 1987.

Brogan, G.X., Friedman, S., McCluskey, C., Cooling, D.S., Berrutti, L., Thode, H.C., and Bock, J.L. Evaluation of a New Quantitative Immunoassay for Serum Myoglobin Versus CK-MB for Ruling Out Acute Myocardial Infarction in the Emergency Department. *Ann. Emerg. Med.* 24: 665-671, 1994.

Juronen, E.I., Viikmaa, M.H. and Mikelsaar, A-V. N. Rapid, Simple and Sensitive Antigen Capture ELISA for the Quantitation of Myoglobin Using Monoclonal Antibodies. *J. Immuno. Met.* 111: 109 - 115, 1988.

Apple, F.S. Acute Myocardial Infarction and Coronary Reprofusion: Serum Cardiac Markers for the 1990s. *Am. J. Clin. Path.* 97: 217-226, 1992.

Mainard, F., Massoubre, B., LeMarec, H., and Madec, Y. Study of a Myoglobin Test in Patients Hospitalized for Suspected Myocardial Infarction. *Clin. Chim. Act.* 153: 1-8, 1985.

Laure, C., Calzolari, C., Bertinchant, J-P, Leclercq, F., Grolleau, R., and Pau, B. Cardiac Specific Immunoenzymometric Assay for Troponin I in the Early Phase of Acute Myocardial Infarction. *Clin. Chem.* 39: 972-979, 1993.

Adams, J.E., Schechtman, K.D., Landt, Y., Ladenson, J.H., and Jaffe, A.S. Comparable Detection of Acute Myocardial Infarction by Creatine Kinase MB Isoenzyme and Cardiac Troponin I. *Clin. Chem.* 40: 1291-1295, 1994.

Adams, J.E., Sicard, G.A., Allen, B.T., Bridwell, K.H., Lenke, L.G., Davila-Roman, V.G., Bodor, G.S., Ladenson, L.H., and Jaffe, A.S. Diagnosis of Perioperative Myocardial Infarction with Measurement of Cardiac Troponin I. *N. Eng. J. Med.* 330: 670-674, 1994.

Brogan, G.X., Hollander, J.E., McCuskey, C.F., Thode, Jr., H.C., Sama, A., Bock, J.L., and the Biochemical Markers for Acute Myocardial Ischemia Study Group. Evaluation of a New Assay for Cardiac Troponin I vs Creatine Kinase-MB for the Diagnosis of Acute Myocardial Infarction. *Acad. Emerg. Med.* 4: 6-12, 1997.

Davis, C.P., Barnett, K., Torre P, and Wacasey, K. Serial Myoglobin Levels for Patients with Possible Myocardial Infarction. *Acad. Emerg. Med.* 3: 590-597, 1996.

Gibler, W.B., Gibler, C.D., Weinshenker, E., Abbottsmith, C., Hedges, J.R., Barsan, W.G., Sperling, M., Chen, I-W., Embry, S., and Kereiakes, D. Myoglobin as an Indicator of Acute Myocardial Infarction. *Ann. Emerg. Med.* 16: 851-856, 1987.

Tucker, J.F., Collins, R.A., Anderson, A.J., Hess, M., Farley, I.M., Hegemann, D.A., Harkins H.J., and Zwicke, D. Value of Serial Myoglobin Levels in the Early Diagnosis of Patients Admitted for Acute Myocardial Infarction. *Ann. Emerg. Med.* 24: 704-708, 1994.

Adams, J.E., Bodor, G., D-Roman, V.G., Delmez, J.A., Apple, F.S., Ladenson J.H., and Jaffe, A.S. Cardiac Troponin I: A Marker with High Specificity for Cardiac Injury. *Circulation* 88: 101-106, 1993.

Buechler, K.F., and McPherson, P.H. Novel Methods for the Assay of Troponin I and T and Complexes of Troponin I and T and Selections of Antibodies for Use in Immunoassays. International Patent WO 96/33415, 18 April, 1995.

Katruckha, A.G., Bereznikova, A.V., Esakova, T.V., Petterson, K., Lövgren, T., Severina, M.E., Pulkki, K., Vuopio-Pulkki, L.-M., and Gusev, N.B. Troponin I is released in bloodstream of patients with acute myocardial infarction not in free form but as complex. *Clin. Chem.* 43: 1379-1385, 1997.

Galvani, M., Ferrini, D., Ghezzi, F., and Ottani, F, Cardiac markers and risk stratification: an integrated approach. *Clin Chim Acta* 311: 9-17, 2001

Meyer, T., Binder, L., Graeber, T., Luthe, H., Kreuzer, H., Oellerich, M., Buchwald, A.B., Superiority of combined CK-MB and troponin I measurements for the early risk stratification of unselected patients presenting with acute chest pain. *Cardiology* 90: 286-294, 1998

de Winter, R.J., Risk stratification with cardiac troponin I in acute coronary syndromes. *J. Am. Coll. Cardiol.* 36: 1824-1826, 2000

Newby, L.K., Storrow, A.B., Gibler, W.B., Garvey, J.L., Tucker, J.F., Kaplan, A.L., Schreiber, D.H., Tuttle, R.H., McNulty, S.E., and Ohman, E.M., Bedside multimarker testing for risk stratification in chest pain units: the chest pain evaluation by creatine kinase-MB, myoglobin, and troponin I (CHECKMATE) study. *Circulation* 103: 1832-1837, 2001.

## DECLARATION OF CONFORMITY

We, Biosite Incorporated,  
11030 Roselle Street  
United States of America

Declare under our sole responsibility that the product:

**Triage® CARDIAC PANEL FOR RAPID  
QUANTIFICATION OF CK-MB,  
MYOGLOBIN, AND TROPONIN I**

To which this declaration relates, is in conformity with the requirements of the  
In-vitro Diagnostic Medical Device Directive (98/79/EC), Annex III.

Signed: *R.H. Dowling*

Title: Director of Quality Assurance and Quality Control

**EC REP** Authorized Representative

Biosite Europe  
Biosite Incorporated Liaison Office  
BP 264  
41/41 BIS Avenue de l'Europe  
78 147 Velizy Cedex FRANCE  
33 134 630665

Intentionally blank

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© 2004, Biosite Incorporated